

**Piezotome 2 Basic Package With 2 Tip Kits**



- Piezotome 2 Console
- 1 Surgical Handpiece
- Any Two Tip Kits
- 1 Large Sterilization Cassette
- 2 Sets Of Autoclavable Tubing
- 20 Irrigation Bag Spikes

PIEZO2-BASIC . . . . . 13,000.00

*Your Choice Of Any Two Tip Kits*

**Piezotome 2 "1 Handpiece" Package With 22 Tips**



- Piezotome 2 Console
- 1 Surgical Handpiece
- Set of 6 Bone Surgery Tips
- Set of 5 Sinus Lift Tips
- Set Of 6 Power Periotome Tips
- Set Of 5 Intralift Tips
- 1 Large Sterilization Cassette
- 2 Sets Of Autoclavable Tubing
- 20 Irrigation Bag Spikes

PIEZO2-1HP . . . . . 14,500.00

**Piezotome 2 "2 Handpiece" Package With 22 Tips**



- Piezotome 2 Console
- 2 Surgical Handpieces
- Set of 6 Bone Surgery Tips
- Set of 5 Sinus Lift Tips
- Set Of 6 Power Periotome Tips
- Set Of 5 Intralift Tips
- 1 Large Sterilization Cassette
- 2 Sets Of Autoclavable Tubing
- 20 Irrigation Bag Spikes

PIEZO2-2HP . . . . . 16,500.00

**SALVIN®**  
Dental Specialties

**NEW!**

**Piezotome 2**  
Piezo Ultrasonic Surgery

**Now With 3X  
More Power!**

**Now With Built In  
LED Lighting!**



**Dr. John Russo**  
Periodontist, Sarasota, FL  
"There is nothing like the Piezotome for lifting the sinus membrane, harvesting autogenous bone, & removing lingual tori while protecting the soft tissue in the floor of the mouth."



**Dr. Michael Sonick**  
Periodontist, Fairfield, CT  
"The use of the Piezotome has been the most revolutionary development in extractions in my career. Most extractions are completed in a matter of minutes with NO loss of alveolar bone."

**Safety**

Cuts bone while preserving soft tissue

**Precision**

Precise cuts to preserve bone



*Piezotome vs. Bur Cut*

**Speed**

60 watts of power for faster cuts  
2 handpiece unit for flexibility during surgery

**Vision**

LED handpiece for excellent visibility



**Everything For Your Implant Practice But The Implants®**  
Salvin Dental Specialties, Inc • Toll Free US & Canada 800-535-6566 • www.salvin.com

© 2009 Salvin Dental Specialties, Inc. All Rights Reserved

11-5-09

© 2009 Salvin Dental Specialties, Inc. All Rights Reserved

**NEW!**


### LED Handpieces

- 100,000 LUX - brighter than sunlight.
- Clear, cool light provides excellent visibility.
- Easy to clean and fully autoclavable at 134°C.



### Selectivity Of The Cut:

The modulated piezo signal (alternation between high and low amplitudes) allows soft tissue relaxation for a clean cut and better healing.



**NEW!**

### Color, Touch-Sensitive Screen

- Backlit, easy to read screen
- Change and save settings at the push of a button



**NEW!**


### PiezoTouch™

The PiezoTouch™ progressive footswitch provides total control of the unit through the footswitch without touching the device screen.

With this new feature, the power of the ultrasonics can be fine tuned where nerves, arteries or membranes are present.



## Bone Surgery




### Ramus Bone Harvesting and Grafting

- BS1 Tip**  
Chin Graft  
#PIEZO2-TIP-BS1
- BS2L Tip**  
Left Ramus Graft  
#PIEZO2-TIP-BS2L
- BS2R Tip**  
Right Ramus Graft  
#PIEZO2-TIP-BS2R
- BS4 Tip**  
Bone Harvesting  
#PIEZO2-TIP-BS4
- BS5 Tip**  
Ridge Expansion  
#PIEZO2-TIP-BS5
- BS6 Tip**  
Osteoplasty  
#PIEZO2-TIP-BS6



## Sinus Lift




### Lateral Sinus Lift

- SL1 Diamond Tip**  
Sinus Window  
#PIEZO2-TIP-SL1
- SL2 Diamond Tip**  
Smooth Edges  
#PIEZO2-TIP-SL2
- SL3 Tip**  
Detach Membrane  
#PIEZO2-TIP-SL3
- SL4 Tip**  
Sinus Elevator  
#PIEZO2-TIP-SL4
- SL5 Tip**  
Sinus Elevator  
#PIEZO2-TIP-SL5



## INTRALIFT™




### Sinus Lift by the Crestal Approach

- TKW1 Tip**  
1.35mm  
#PIEZO2-TIP-TKW1
- TKW2 Tip**  
2.10mm  
#PIEZO2-TIP-TKW2
- TKW3 Tip**  
2.35mm  
#PIEZO2-TIP-TKW3
- TKW4 Tip**  
2.80mm  
#PIEZO2-TIP-TKW4
- TKW5 Tip**  
Irrigated Tip  
#PIEZO2-TIP-TKW-5



## EXTRACTION



### Extraction and Immediate Placement

- LC1-90 Tip**  
90° Perioste Tip  
#PIEZO2-TIP-LC1-90
- LC1 Tip**  
Perioste Tip  
#PIEZO2-TIP-LC1
- LC2 Tip**  
Perioste Slim Tip  
#PIEZO2-TIP-LC2
- LC2L Tip**  
45° Left Slim Tip  
#PIEZO2-TIP-LC2L
- LC2R Tip**  
45° Right Slim Tip  
#PIEZO2-TIP-LC2R
- Ninja Tip**  
Tooth Sectioning Tip  
#PIEZO2-TIP-NINJA



## CROWN EXTENSION



### NEW! Crown Extension

- CE1 Tip**  
Osteoplasty  
#PIEZO2-TIP-CE1
- CE2 Tip**  
Interproximal  
#PIEZO2-TIP-CE2
- CE3 Tip**  
Around Root  
#PIEZO2-TIP-CE3
- CE4 Tip**  
Bone Reshaping  
#PIEZO2-TIP-CE4

